


# How to Connect Your Android Phone to Your Widex EVOKE Hearing Aids

## Step by step guide

**NOTE:** Red circles indicate where to 'touch/click on your Phone'

**NOTE** Android Phone require an intermediary device called the [Widex COM-DEX](#) to stream phone calls and media (music, podcast, videos) from your phone to your hearing aids.


**1.** Ensure **Bluetooth** on your phone is activated.



2. Click on the **Settings App** on your Android phone



3. Under the **Network** subfolder, click on **Bluetooth**




**4. Now, turn your hearing aids off and then on** by opening and closing the battery door.


**5. Wait 1-2 minutes** whilst the Phone automatically 'searches' for your hearing aids.

**6. Once found, your hearing aids will be shown in the list under Available Devices.** In our example, our hearing aids are titled 'test Hearing aids'

Select each hearing device to **pair** to your Phone.




**6. They should now appear under the list Paired Devices**



7. Once paired, exit the Bluetooth menu by pressing the **home button**

8. Find your **Widex EVOKE App** and click it



9. Congratulations, your hearing aids are now successfully paired to your Phone. You can now:

- change your hearing aid programs
- adjust multiple parameters of your hearing aids through the convenience of your Phone
- receive phone calls (requires Com-Dex\*)
- stream music, podcasts and video sound (requires Com-Dex\*)

\*Android Phone require an intermediary device called the [Widex COM-DEX](#) for these features to work between your Android phone and your hearing aids.

